

05/03/2020

3M N95 Masks Treatable by Battelle CCDS™



3M™ 1860(S)



3M1860



3M™ 8000 & 8205



3M™ Aura 1870+



3M™ Aura 9210+/37192



3M™ Aura 9211+/37193



3M™ 8110S & 8210
& 8210plus



3M™ 8210V



3M™ 8211



3M™ 8511



3M™
8515/07189



3M™ 8516

Other N95 Masks Treatable by Battelle CCDS™



Honeywell N1105



Honeywell N1125



Honeywell NBW95




Kimberly-Clark 46727,
46767, and 46867



Kimberly-Clark
62126 and 62355



Moldex 1500 Series



Moldex 2200



Gerson 2140C



Gerson 1730

Contain Cellulose – **DO NOT TREAT**



Check for Cellulose Acetate in undefined masks by placing filter, shell, and coverweb (not straps or nose foam) in Butyric Acid.

Polypropylene and Polyester have satisfactory chemical resistance to it, but Cellulose Acetate (CA) does not so degradation should be evident if the mask contains CA.